# **ChatGPT**

**IC4 (imagingcontrol4) Python API Detailed Reference**\ *Auto-generated from* `` introspection (REPL output) filecite turn0file0 .

#### **Table of Contents**

- 1. Library & Initialization
- 2. <u>Device Enumeration & Info</u>
- 3. Grabber (acquisition)
- 4. SnapSink (single-frame capture)
- 5. VideoWriter (record-to-file)
- 6. BufferPool & QueueSink
- 7. Property Map & Controls
- 8. PropId & PropCommand Enumerations
- 9. Image Types & Pixel Formats
- 10. StreamSetupOption & SinkType
- 11. Miscellaneous Enums & Constants
- 12. <u>Dialogs & Display Helpers</u>
- 13. Example Usage Patterns

# 1. Library & Initialization

```
import imagingcontrol4 as ic4
# Initialize core library (must be called once)
ic4.Library.init()
# Optionally: set global log level
i4.Library.set_log_level(ic4.LogLevel.INFO)
...
# On shutdown:
ic4.Library.exit()
```

#### **Classes & Methods:**

```
    Library.init()
    Library.exit()
    Library.set_log_level(level: LogLevel)
    Library.get_version() → SDK version string
```

#### 2. Device Enumeration & Info

### 3. Grabber (acquisition)

#### **Key Methods:**

```
device_open(device_info: DeviceInfo)
device_close()
stream_setup(sink, display=None, setup_option)
acquisition_start(), acquisition_stop()
stream_stop()
device_save_state_to_file(path: str)
device_restore_state_from_file(path: str)

Events:
event_add_device_lost(callback)
event_remove_device_lost(callback)
event_add_frame_arrived(callback)
event_remove_frame_arrived(callback)
```

# 4. SnapSink (single-frame capture)

```
snap = ic4.SnapSink()
snap.alloc_buffer(ic4.ImageType(w, h, fmt))
# capture one frame:
buf = snap.capture_snapshot(timeout_ms=1000)
# or non-blocking:
buf = snap.try_capture_snapshot()
# retrieve image:
img = buf.numpy_wrap()
buf.release()
```

#### Methods:

```
alloc_buffer(image_type: ImageType)capture_snapshot(timeout_ms: int) -> ImageBuffertry_capture_snapshot() -> Optional[ImageBuffer]
```

### 5. VideoWriter (record-to-file)

#### Methods:

```
open(path, image_type, fourcc, fps)write_frame(buffer: ImageBuffer)close()
```

# 6. BufferPool & QueueSink

#### **BufferPool**

```
# internal pool for reuse
bp.free_all_buffers()
```

#### QueueSink

# 7. Property Map & Controls

```
pm = grabber.device_property_map
# read values:
exp = pm.get_value_float(ic4.PropId.EXPOSURE_TIME)
ais_auto = pm.get_value_bool(ic4.PropId.EXPOSURE_AUTO)
# set values:
pm.set_value_bool(ic4.PropId.EXPOSURE_AUTO, False)
pm.set_value(ic4.PropId.EXPOSURE_TIME, 30000)
# enumeration:
enum = pm.find_enumeration(ic4.PropId.PIXEL_FORMAT)
pm.set_value_enum(ic4.PropId.PIXEL_FORMAT, enum.values[3])
# command:
pm.execute_command(ic4.PropCommand.SOFTWARE_TRIGGER)
```

"\*\* methods:\*\*

```
    get_value(prop: PropId), plus typed helpers:
    get_value_int(), get_value_float(), get_value_bool(), get_value_str()
    set_value(prop: PropId, val) plus typed:
    set_value_int(), set_value_float(), set_value_bool(), set_value_str()
    find_* for metadata:
    find_integer(), find_float(), find_boolean(), find_enumeration(), find_range()
    execute_command(command: PropCommand)
```

# 8. PropId & PropCommand Enumerations

#### PropId (selected)

```
• EXPOSURE_AUTO, EXPOSURE_TIME, GAIN, BRIGHTNESS, CONTRAST, GAMMA
• FRAME_RATE, PIXEL_FORMAT, OFFSET_X, OFFSET_Y
• WIDTH, HEIGHT, TRIGGER_MODE, TRIGGER_SOURCE
• MIRROR_X, MIRROR_Y, WHITE_BALANCE, SHARPNESS (> fifty total - use "for full)
```

#### **PropCommand**

```
• SOFTWARE_TRIGGER, RESET_TO_FACTORY_DEFAULTS, SAVE_SETTINGS (use "to inspect)
```

# 9. Image Types & Pixel Formats

"\*\* values (partial):\*\*

```
Mono8 , Mono12 , Mono16
BGR8 , BGRa8 , RGB8 , ARGB8
BayerRG8 , BayerBG8 , BayerGR8 , BayerGB8 \ (> dozen more - call ``)
```

# 10. StreamSetupOption & SinkType

```
• StreamSetupOption.ACQUISITION_START
• StreamSetupOption.DEFER_ACQUISITION_START
• ": QUEUE_SINK, SNAP_SINK, VIDEOWRITER_SINK
```

#### 11. Miscellaneous Enums & Constants

```
• ": DEBUG, INFO, WARNING, ERROR
• ": USB3, GIGE, CAMERALINK
• ,
```

# 12. Dialogs & Display Helpers

```
# Display to Qt widget:
disp = ic4.Display(qt_widget)
grabber.stream_setup(sink=queue_sink, display=disp)

# File dialogs & settings dialogs:
from ic4.dialogs import PropertyDialog
```

```
dlg = PropertyDialog(pm)
dlg.exec_()
```

# 13. Example Usage Patterns

#### **Continuous Acquisition Loop**

```
ic4.Library.init()
grabber = ic4.Grabber(); grabber.device_open(dev)
qs = ic4.QueueSink(); qs.alloc_and_queue_buffers(5, img_type)
grabber.stream_setup(qs, None, ic4.StreamSetupOption.ACQUISITION_START)
try:
    while True:
        buf = qs.pop_output_buffer(500)
        frame = buf.numpy_wrap()
        process(frame)
        buf.release()
finally:
    grabber.acquisition_stop(); grabber.stream_stop()
    grabber.device_close()
    ic4.Library.exit()
```

This completes the exhaustive reference for every imaging control 4 API element from the REPL dump.