

## # PRIMAcquisition

**\*\*PRIMAcquisition\*\*** (PRIM) is a Python-based application for synchronized acquisition of pressure data (from an Arduino-controlled pressure transducer) and live camera imaging. Designed for vascular physiology experiments, PRIMAcquisition displays live pressure traces, previews a high-speed camera feed, and saves synchronized recordings (as a CSV and TIFF stack) for offline analysis.

---

## ## 🚀 Quick Start – Run an Experiment

- 1 **\*\*Connect PRIM Device\*\*** – Select the Arduino COM port and click **\*\*Connect PRIM Device\*\***
- 2 **\*\*Set Up Camera\*\*** – Choose camera & resolution → click **\*\*Start Camera\*\***
- 3 **\*\*Adjust Exposure/Gain\*\*** – Use controls to fine-tune camera settings
- 4 **\*\*Zero PRIM\*\*** – Make sure pressure is at zero
- 5 **\*\*Start Recording\*\*** – Click **\*\*Start Recording\*\*** to begin acquisition
- 6 **\*\*Stop Recording\*\*** – Click **\*\*Stop Recording\*\*** when finished
- 7 **\*\*Playback & Export\*\*** – Click **\*\*Playback\*\*** to review, export a frame (PNG) or full TIFF with overlay

## ## Features

### ### 🎥 Real-Time Pressure + Video Recording

- **\*\*Arduino is the master clock\*\***:
  - Generates precise trigger pulses (``CamTrig``) to the camera for each frame.
  - Immediately sends synchronized serial data: ``frame_index, elapsed_time_s, pressure_value``.
- **\*\*App listens for first Arduino tick\*\***, then begins writing data.
- **\*\*Perfect sync\*\*** is maintained:
  - Each TIFF frame corresponds to exactly one Arduino trigger.
  - Each pressure value corresponds to exactly one recorded frame.
  - The matching pressure value is stored in each TIFF frame's metadata.

### ### 📈 Live Pressure Plotting

- Receives pressure readings from Arduino (115200 baud).
- Displays live pressure trace in real time.
- Shows frame index, elapsed time, and pressure in the top panel.

### ### 📷 High-Speed Camera Preview & Control

- Integrates with The Imaging Source DMK cameras via IC Imaging Control 4 (IC4).
- Lists connected USB3 Vision cameras with supported resolutions.
- Displays a full-resolution preview using OpenGL viewfinder.
- Camera settings (exposure, gain, brightness) adjustable via sliders.

### ### 📁 Synchronized Recording Output

- Folder structure:

PRIM\_ROOT/YYYY-MM-DD/FillN/

```
| recording_*.csv      ← Pressure + frame index  
| recording_*.tif      ← Grayscale video, one frame per Arduino trigger  
|`
```

- Files are automatically named and saved in time-stamped subfolders.
- Environment variable `PRIM\_RESULTS\_DIR` overrides default save location.
- After recording you can open **\*\*Playback\*\*** to view the TIFF with pressure overlay (toggleable), export an annotated copy, or save a single-frame snapshot.

----

### ## 🛠 Under the Hood

PRIMAcquisition uses a **\*\*hardware-triggered acquisition model\*\***:

| Component          | Role                                                      |
|--------------------|-----------------------------------------------------------|
| <b>**Arduino**</b> | 🕒 Master clock; sends trigger pulses and serial data      |
| <b>**Camera**</b>  | 📷 Triggered via `CamTrig` pin from Arduino                |
| <b>**App**</b>     | 🧠 Waits for first Arduino message, then saves video + CSV |

Each cycle:

1. Arduino triggers a camera frame via digital HIGH → LOW on `CamTrig`.
2. Arduino immediately sends pressure data and timestamp via serial.
3. App receives both the frame and serial line, saves them together.

Result: **\*\*pixel-perfect and time-accurate alignment\*\*** of pressure + video frames.

- **\*\*SerialThread\*\*** - Reads Arduino serial data and emits `(frameIndex, timestamp, pressure)`.
- **\*\*SDKCameraThread\*\*** - Opens IC4 camera, configures settings, streams frames to the GUI.
- **\*\*RecordingManager\*\*** - Writes synchronized CSV and TIFF files in a

background thread.

– **PlaybackWindow** – Reloads recorded files and overlays pressure data on frames.

Recording starts when the first Arduino tick is received. Each pressure value is embedded in its corresponding TIFF frame metadata, ensuring **perfect synchronization**.

---

## ## Installation

### 1. **Clone the Repository**

```
```bash
git clone https://github.com/your-repo/PRIMAcquisition.git
cd PRIMAcquisition/prim_app
```
```

### 2. **Create a Virtual Environment**

```
```bash
python -m venv .venv
.venv\Scripts\activate    # On Windows
```
```

### 3. **Install Dependencies**

```
```bash
pip install -r requirements.txt
```
```

Or manually:

```
```bash
pip install pyqt5 imagingcontrol4 pyserial numpy tifffile
```
```

### 4. **Install IC4 SDK**

- \* Required for DMK camera support.
- \* Ensure GenTL Producer is installed and camera enumerates via IC4.

---

## ## Running the App

### 1. **Start App**

```
```bash
python prim_app.py
```
```

```

## 2. \*\*Connect Arduino\*\*

- \* Select COM port (e.g., COM8), then click \*\*Connect PRIM Device\*\*.
- \* Confirm status shows "Connected".

## 3. \*\*Configure Camera\*\*

- \* Select camera and resolution.
- \* Click \*\*Start Camera\*\* to preview live feed.

## 4. \*\*Start Recording\*\*

- \* Use menu: \*\*Acquisition → Start Recording\*\* or press \*\*Ctrl+R\*\*.
- \* App waits for first Arduino `T...` tick to begin synchronized recording.
- \* TIFF and CSV are saved in the correct folder.

## 5. \*\*Stop Recording\*\*

- \* Use menu: \*\*Acquisition → Stop Recording\*\* or press \*\*Ctrl+T\*\*.
- \* Files are finalized and closed cleanly.

---

## ## Arduino Firmware (PRIM\_v3\_02)

- \* Loops every `startup.timeDelay` milliseconds.
- \* Sends a `CamTrig` HIGH-LOW pulse to trigger one camera frame.
- \* Immediately sends formatted serial data:

```

```
frame_index, elapsed_time_s, pressure_value
```

```

- \* Operates with precise timing via `micros()` and `millis()`.

---

## ## Troubleshooting

| Issue                | Fix                                         |
|----------------------|---------------------------------------------|
| Camera not listed    | Verify IC4 SDK + GenTL Producer installed   |
| Serial shows no data | Check Arduino COM port and baud rate        |
| TIFF won't open      | Use ImageJ/Fiji or Python `tifffile`        |
| Dropped frames       | Use USB 3.0 and reduce resolution if needed |

---

## ## License

Creative Commons Attribution-NonCommercial-ShareAlike 4.0  
International

[<https://creativecommons.org/licenses/by-nc-sa/4.0/>](<https://creativecommons.org/licenses/by-nc-sa/4.0/>)